

Product Datasheet

Mouse CRH ELISA Kit (orb438842)

Description	Mouse Corticotropin Releasing Hormone (CRH) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8CIT0
Sensitivity	4.45pg/mL
Entrez	12918
Expiration Date	Please enquire.