

Product Datasheet

Rat CRHR1 ELISA Kit (orb438844)

Description	Rat Corticotropin Releasing Hormone Receptor 1 (CRHR1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35353
Sensitivity	0.123ng/mL
Entrez	58959
Expiration Date	Please enquire.