

Product Datasheet

Rat CRHR2 ELISA Kit (orb438846)

Description	Rat Corticotropin Releasing Hormone Receptor 2 (CRHR2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P47866
Sensitivity	0.134ng/mL
Entrez	64680
Expiration Date	Please enquire.