

Product Datasheet

Rat CRH ELISA Kit (orb438847)

Description	Rat Corticotropin Releasing Hormone (CRH) competitive inhibition ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01143
Sensitivity	4.58pg/mL
Entrez	81648
Expiration Date	Please enquire.