

Product Datasheet

Human CRIM1 ELISA Kit (orb438849)

Description	Human Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NZV1
Sensitivity	0.133ng/mL
Entrez	51232
Expiration Date	Please enquire.