

Product Datasheet

Human CRISP3 ELISA Kit (orb438852)

Description	Human Cysteine Rich Secretory Protein 3 (CRISP3) Sandwich ELISA Kit
Reactivity	Human
Range	1.875-120ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva, sweat and other biological fluids.
Sample Types	serum, plasma, saliva, sweat and other biological fluids.
Uniprot ID	P54108
Sensitivity	0.68ng/mL
Entrez	10321
Expiration Date	Please enquire.