

Product Datasheet

Mouse CRLF1 ELISA Kit (orb438856)

Description	Mouse Cytokine Receptor Like Factor 1 (CRLF1) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9JM58
Sensitivity	5.3pg/mL
Entrez	12931
Expiration Date	Please enquire.