

Product Datasheet

Chicken CRP ELISA Kit (orb438863)

Description	Chicken C Reactive Protein (CRP) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2EJU6
Sensitivity	0.211ng/mL
Entrez	429786
Expiration Date	Please enquire.