

Product Datasheet

Guinea pig CRP ELISA Kit (orb438864)

Description	Guinea pig C Reactive Protein (CRP) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P49254
Sensitivity	7.4pg/mL
Entrez	100727100
Expiration Date	Please enquire.