

Product Datasheet

Human CRTC3 ELISA Kit (orb438872)

Description	Human CREB Regulated Transcription Coactivator 3 (CRTC3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6UUV7
Sensitivity	0.112ng/mL
Entrez	64784
Expiration Date	Please enquire.