

Product Datasheet

Human CRYaB ELISA Kit (orb438878)

Description	Human Crystallin Alpha B (CRYaB) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P02511
Sensitivity	0.216ng/mL
Entrez	1410
Expiration Date	Please enquire.