

Product Datasheet

Human CSN1 ELISA Kit (orb438884)

Description	Human Casein Alpha (CSN1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates, breast milk and other biological fluids.
Sample Types	tissue homogenates, cell lysates, breast milk, cell culture supernatants and other biological fluids.
Uniprot ID	P47710
Sensitivity	0.55ng/mL
Entrez	1446
Expiration Date	Please enquire.