

Product Datasheet

Bovine CSN2 ELISA Kit (orb438885)

Description	Bovine Casein Beta (CSN2) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	milk.
Sample Types	milk.
Uniprot ID	P02666
Sensitivity	27ng/mL
Entrez	281099
Expiration Date	Please enquire.