

Product Datasheet

Human CSN2 ELISA Kit (orb438886)

Description	Human Casein Beta (CSN2) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	breast milk and other biological fluids.
Sample Types	breast milk and other biological fluids.
Uniprot ID	P05814
Sensitivity	27ng/mL
Entrez	1447
Expiration Date	Please enquire.