

Product Datasheet

Bovine CSN3 ELISA Kit (orb438887)

Description	Bovine Casein Kappa (CSN3) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.312-20µg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	milk.
Sample Types	milk.
Uniprot ID	P02668
Sensitivity	0.101µg/mL
Entrez	281728
Expiration Date	Please enquire.