



Product Datasheet

CD41 antibody (orb43889)


Description	Mouse monoclonal antibody to CD41
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP
Immunogen	Bovine trombocytes.
Target	CD41/CD61
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD41/CD61 positive b...



Separation of bovine CD41/CD61 positive ...



Flow cytometry surface staining pattern ...