

Product Datasheet

Rat CST3 ELISA Kit (orb438896)

Description	Rat Cystatin 3 (CST3) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	P14841
Sensitivity	0.068ng/mL
Entrez	25307
Expiration Date	Please enquire.