

Product Datasheet

Human CST6 ELISA Kit (orb438898)

Description	Human Cystatin 6 (CST6) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, sweat, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, sweat, cell culture supernatants and other biological fluids.
Uniprot ID	Q15828
Sensitivity	0.053ng/mL
Entrez	1474
Expiration Date	Please enquire.