

Product Datasheet

Mouse CSTA ELISA Kit (orb438900)

Description	Mouse Cystatin A (CSTA) Sandwich ELISA Kit
Reactivity	Mouse
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P56567
Sensitivity	14.25pg/mL
Entrez	209294
Expiration Date	Please enquire.