

Product Datasheet

Rat CT1 ELISA Kit (orb438904)

Description	Rat Cardiotrophin 1 (CT1) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q63086
Sensitivity	5.3pg/mL
Entrez	29201
Expiration Date	Please enquire.