

Product Datasheet

Mouse CTSK ELISA Kit (orb438934)

Description	Mouse Cathepsin K (CTSK) Sandwich ELISA Kit
Reactivity	Mouse
Range	12.5-800pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P55097
Sensitivity	5.4pg/mL
Entrez	13038
Expiration Date	Please enquire.