

Product Datasheet

Rat CTSK ELISA Kit (orb438936)

Description	Rat Cathepsin K (CTSK) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	035186
Sensitivity	6.5pg/mL
Entrez	29175
Expiration Date	Please enquire.