

Product Datasheet

Mouse CTSL ELISA Kit (orb438940)

Description	Mouse Cathepsin L (CTSL) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q543M3
Sensitivity	23.7pg/mL
Entrez	13039
Expiration Date	Please enquire.