

Product Datasheet

Porcine CTSL ELISA Kit (orb438941)

Description	Porcine Cathepsin L (CTSL) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	E1UJ48
Sensitivity	0.132ng/mL
Entrez	100515919
Expiration Date	Please enquire.