

Product Datasheet

Human CTSV ELISA Kit (orb438946)

Description	Human Cathepsin V (CTSV) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O60911
Sensitivity	10.9pg/mL
Entrez	1515
Expiration Date	Please enquire.