

Product Datasheet

Canine CX43 ELISA Kit (orb438969)

Description	Canine Connexin 43 (CX43) Sandwich ELISA Kit
Reactivity	Canine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	J9P7X9
Sensitivity	0.074ng/mL
Entrez	403418
Expiration Date	Please enquire.