

Product Datasheet

Human CXCL16 ELISA Kit (orb438976)

Description	Human Chemokine C-X-C-Motif Ligand 16 (CXCL16) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9H2A7
Sensitivity	0.067ng/mL
Entrez	58191
Expiration Date	Please enquire.