

Product Datasheet

Rat CXCR3 ELISA Kit (orb438980)

Description	Rat Chemokine C-X-C-Motif Receptor 3 (CXCR3) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9JII9
Sensitivity	0.059ng/mL
Entrez	84475
Expiration Date	Please enquire.