

Product Datasheet

Human CXCR7 ELISA Kit (orb438986)

Description	Human Chemokine C-X-C-Motif Receptor 7 (CXCR7) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P25106
Sensitivity	0.062ng/mL
Entrez	57007
Expiration Date	Please enquire.