

Product Datasheet

Mouse CYP11B1 ELISA Kit (orb438995)

Description	Mouse Cytochrome P450 11B1 (CYP11B1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q3TG86
Sensitivity	0.064ng/mL
Entrez	110115
Expiration Date	Please enquire.