

Product Datasheet

Rat CYP1A1 ELISA Kit (orb438998)

Description	Rat Cytochrome P450 1A1 (CYP1A1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P00185
Sensitivity	0.129ng/mL
Entrez	24296
Expiration Date	Please enquire.