

Product Datasheet

Monkey CYP1A1 ELISA Kit (orb438999)

Description	Monkey Cytochrome P450 1A1 (CYP1A1) Sandwich ELISA Kit
Reactivity	Monkey
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q6GUR1
Sensitivity	0.056ng/mL
Entrez	678694
Expiration Date	Please enquire.