

Product Datasheet

Human CYP1A2 ELISA Kit (orb439000)

Description	Human Cytochrome P450 1A2 (CYP1A2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05177
Sensitivity	0.051ng/mL
Entrez	1544
Expiration Date	Please enquire.