

Product Datasheet

Human CYP27A1 ELISA Kit (orb439005)

Description	Human Cytochrome P450 27A1 (CYP27A1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q02318
Sensitivity	0.061ng/mL
Entrez	1593
Expiration Date	Please enquire.