

Product Datasheet

Human CYP2C18 ELISA Kit (orb439009)

Description	Human Cytochrome P450 2C18 (CYP2C18) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P33260
Sensitivity	0.247ng/mL
Entrez	1562
Expiration Date	Please enquire.