

Product Datasheet

Human CYP2E1 ELISA Kit (orb439014)

Description	Human Cytochrome P450 2E1 (CYP2E1) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05181
Sensitivity	0.29ng/mL
Entrez	1571
Expiration Date	Please enquire.