

Product Datasheet

Mouse CYP7A1 ELISA Kit (orb439020)

Description	Mouse Cytochrome P450 7A1 (CYP7A1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q64505
Sensitivity	0.139ng/mL
Entrez	13122
Expiration Date	Please enquire.