

Product Datasheet

Human CYR61 ELISA Kit (orb439028)

Description	Human Cysteine Rich Protein, Angiogenic Inducer 61 (CYR61) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	000622
Sensitivity	13.4pg/mL
Entrez	3491
Expiration Date	Please enquire.