

Product Datasheet

Human DAG1 ELISA Kit (orb439040)

Description	Human Dystrophin Associated Glycoprotein 1 (DAG1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q14118
Sensitivity	0.086ng/mL
Entrez	1605
Expiration Date	Please enquire.