

Product Datasheet

Rat DAO ELISA Kit (orb439046)

Description	Rat Diamine Oxidase (DAO) Sandwich ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O35078
Sensitivity	1.33ng/mL
Entrez	114027
Expiration Date	Please enquire.