

Product Datasheet

Mouse DCLK1 ELISA Kit (orb439064)

Description	Mouse Doublecortin Like Kinase 1 (DCLK1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JLM8
Sensitivity	0.059ng/mL
Entrez	13175
Expiration Date	Please enquire.