

Product Datasheet

Human DCN ELISA Kit (orb439065)

Description	Human Decorin (DCN) Sandwich ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07585
Sensitivity	35pg/mL
Entrez	1634
Expiration Date	Please enquire.