

Product Datasheet

Human DCX ELISA Kit (orb439067)

Description	Human Doublecortin (DCX) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O43602
Sensitivity	0.058ng/mL
Entrez	1641
Expiration Date	Please enquire.