

Product Datasheet

Human DDIT3 ELISA Kit (orb439075)

Description	Human DNA Damage Inducible Transcript 3 (DDIT3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P35638
Sensitivity	0.138ng/mL
Entrez	1649
Expiration Date	Please enquire.