

Product Datasheet

Human DEFa3 ELISA Kit (orb439084)

Description	Human Defensin Alpha 3, Neutrophil Specific (DEFa3) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P59666
Sensitivity	16.47pg/mL
Entrez	1668
Expiration Date	Please enquire.