

Product Datasheet

Human DEFb124 ELISA Kit (orb439092)

Description	Human Defensin Beta 124 (DEFb124) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8NES8
Sensitivity	5.6pg/mL
Entrez	245937
Expiration Date	Please enquire.