

Product Datasheet

Human DEFb2 ELISA Kit (orb439103)

Description	Human Defensin Beta 2 (DEFb2) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	015263
Sensitivity	13.55pg/mL
Entrez	1673
Expiration Date	Please enquire.