

Product Datasheet

Mouse DIO3 ELISA Kit (orb439130)

Description	Mouse Deiodinase, Iodothyronine, Type III (DIO3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q91ZI8
Sensitivity	0.115ng/mL
Entrez	107585
Expiration Date	Please enquire.