

Product Datasheet

Human DKK2 ELISA Kit (orb439135)

Description	Human Dickkopf Related Protein 2 (DKK2) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9UBU2
Sensitivity	34pg/mL
Entrez	27123
Expiration Date	Please enquire.