


Product Datasheet

LST1 antibody (orb43914)

Description	Mouse Monoclonal to LST1.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	KLH-fused peptide corresponding to amino acids SSEGPDLRGRDKRGT of human LST1
Target	LST1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunoprecipitation: Positive control: RAJI human Burkitt lymphoma cell line. Western blotting: The antigen migrates on SDS PAGE gels as approximately 25-28 kDa molecule.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	000453
Entrez	7940
Dilution Range	Immunoprecipitation: Positive control: RAJI human Burkitt lymphoma cell line. Western blotting: The antigen migrates on SDS PAGE gels as approximately 25-28 kDa molecule.
Expiration Date	12 months from date of receipt.



Western blotting analysis (non-reducing ...)